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EN User Instruction DC Photoionisation Lamp

Application area and intended use


The DC photoionisation lamp is a gas discharge lamp and has been designed to be used within photoionisation detectors, gas chromatography, mass spectrometry and ion mobility spectrometry devices.

 This DC Photoionisation Lamp is exempted from Regulation (EU) 2019/2020 as amended in accordance with Annex III, 2d and not intended for general lighting. Any use in other than the intended applications is not admissible.

Safety information

 → Follow the safety instructions. Non-observance of the safety notes may cause hazards to persons, the environment and the product.

DANGER dangerous voltage

 Electric shock may occur when working under voltage. Danger to life.

→ Disconnect product from the power supply. Check the system is de-energised. Secure against restarting.

WARNING ultraviolet radiation

 This product belongs to Risk group 1. Low risk, pursuant to EN 62471 (Photobiological safety of lamps and lamp systems). Safe for most applications apart from direct long exposure very close proximity (<1m) viewing. Minimise exposure to eyes and skin. Use appropriate shielding.

In normal applications the lamp is inside an instrument and cannot be viewed directly. Even if the lamp were to be somehow run in the open the optical risks are extremely low. See note above.

Transport and handling

 → Transport the product in its original packing to the place of installation. The cleanliness of the PID lamp crystal window is extremely important to ensure high photon output and good detector sensitivity. Finger grease on the crystal window leads to primary VUV line output losses. Avoid touching the lamp with bare hands. Wear gloves.

→ Prior to installation, check the DC photoionisation lamp for mechanical damage. **Never use damaged DC photoionisation lamps.**

Installation

De-energise the unit before commencing work on it in accordance with the relevant safety precautions for working with electrical equipment. Install or replace the lamp in the instrument according to the instruction manual of the instrument in which the lamp is used. Allow the lamp to cool sufficiently to prevent burns from hot surfaces in the case of heated instrumentation. The lamp itself does not generate much heat as the input power is normally less than 1W, so in non-

heated equipment there are no risks of heat burns. Only use the tools described in the maintenance manual to replace the lamp.

The lamp has bare electrode surfaces and the equipment in which the lamp is designed to fit has the appropriate mounts and electrical contacts.

Disposal

RoHS: Photoionisation lamps are in scope of the Directive 2011/65/EU (RoHS II, Category 5 Lighting equipment) as amended and are RoHS II compliant.

China RoHS: Photoionisation lamps are subject to the People's Republic of China law The Management Methods for the Restriction of the Use of Hazardous Substances in Electrical and Electronic Products, order No. 32 of January 6th, 2016. DC photoionisation lamps must be labelled with the pictured, orange EFUP 50 as stipulated in the Standard SJ/T 11364-2014. The names of the hazardous substances and where in the lamp they are present above the concentration limits are specified in the following table. DC photoionisation lamps should be recycled in accordance with the relevant local or national legal requirements for electrical and electronic waste.

Component Name	Hazardous Substances					
	Lead (Pb)	Mercury (Hg) and Compounds	Cadmium (Cd)	Chromium VI Compounds (Cr(VI))	Polybrominated Biphenyls (PBB)	Polybrominated Diphenyl Ethers (PBDE)
Glass frit	X	0	0	0	0	0

REACH: Note pursuant to Article 33 REACH 1907/2006: the glass frit contains more than 0,1% lead (CAS No. 7439-92-1).

WEEE: DC photoionisation lamps are in scope of Directive 2012/19/EU WEEE and can be classified according to Annex 3(3) and 4(3) (Lamps: Gas Discharge Lamps). DC photoionisation lamps may not be disposed of in the normal commercial or industrial waste. EU users must dispose of the lamps at a communal recycling center for electrical and electronic waste. Other users should recycle, otherwise dispose of the lamps in accordance with the relevant legal requirements.

EU Declaration of Conformity

Heraeus Noblelight declares that this product is in conformity with the following Directives: RoHS 2011/65/EU. Low Voltage 2014/35/EU. If required, a complete Declaration of Conformity can be requested.

Complaints

→ In the event of a complaint, return the DC photoionization lamp in its original packaging to the dealer or the address provided.

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